



Instruments for Composites



A-DAPT™

Composite Instrument

The A-DAPT™ composite instruments facilitate the adaptation, modelling and contouring of composite and glass-ionomer restorations.

- Allow precise and fine contouring and modelling of composite restorations without the inconvenience of the material sticking to the instrument.
- Handles are lightweight and thin as pencils to secure the best possible tactile feel.
- Tips are made from a non-stick silicone material of the perfect softness.
- Handle and tip can be autoclaved by max. 135°C/275°F.

Packing:

- 1 straight stainless steel handle with 2 brush shaped tips.
- 1 angled (60°) stainless steel handle with 2 conical, rounded tips.

Available as single packs or in a set.

Tips are available in sets of two.



RØNVIG
Dental Mfg. A/S

Gl. Vejlevej 59 · DK-8721 Daugaard · Denmark
Tel.: (+45) 70 23 34 11 · Fax: (+45) 76 74 07 98
e-mail: export@ronvig.com · www.ronvig.com