

Instruments for Restoration



SEE-IT™

SEE-IT™ caries indicator and root canal finder is designed to identify and distinguish between carious dentine infused with bacteria and carious dentine which is not infected.

SEE-IT™ is a valuable help in ensuring that all infected dentine is removed while preserving as much as possible of the underlying non-infected (remineralizable) dentine.

- Stains infected carious dentine also in difficult to see locations (e.g. under overhanging enamel) and detects difficult-to-find root canals
- Is an easy-to-handle thixotropic gel (no dripping).
- Is easily applied in the cavity by using a small stick, sponge or mini-brush
- Available in red and blue (10 ml)
- Available in bottles of 10 ml. and in syringes of 1,4 ml. (set of 5 pcs.)



RØNVIG
Dental Mfg. A/S

Gl. Vejlevej 59 · DK-8721 Daugaard · Denmark
Tel.: (+45) 70 23 34 11 · Fax: (+45) 76 74 07 98
e-mail: export@ronvig.com · www.ronvig.com